



KA-1000
Autorefractor

KW-2000

Combined autorefractor and keratometer

Quick and simple for accurate refraction

Kowa's single footprint autorefractors offer quick and easy objective measurements of the refractive error of the eye and provide the data you need to assist with both spectacle and contact lens fitting. The auto start function makes operation straightforward – just focus and measurement begins automatically. Interpretation of the results is simple with an on-screen readout display using coloured icons as well as a paper printout. Small pupils down to $\varnothing 2.3\text{mm}$ can be measured and IOL implants considered. In addition, the KW-2000 has the added function of keratometry with central measurement as well as the periphery in four directions.



Special features

KA-1000 and KW-2000

- Easy auto start operation
- Accurate results
- Colour display
- Paper printout
- Connectivity with third party electronic health records (EHR)

KW-2000

- Central keratometry readings
- Peripheral K-readings in 4 directions
- Refraction and central K-readings in one measurement



On-screen display with
coloured icons



Paper printout